

PC-based portable

# Piccolo



300 grams



USB-powered



- AC, BC and Speech, QuickSIN™
- Weighs less than 300 grams
- Cost-effective solution
- Type 3 pure tone audiometer



PC-based portable  
**Piccolo**



**Piccolo** is a portable computer-controlled audiometer, offering air, bone and speech audiometry functions. A lightweight, compact design, powerful, user-friendly software, and Noah connectivity combine to make Piccolo the ideal choice for professionals on the move.

Available in three versions:

- Speech: pure tone and speech audiometry - AC and BC outputs. Only the Speech model includes a 2 channels Master Hearing Aid and Qui- ckSIN™.
- Basic: pure tone audiometry - only AC output
- Plus: pure tone audiometry - AC and BC outputs



**INCREDIBLY PORTABLE**

Piccolo is extremely light (less than 300 grams), boasts a very small footprint, and does not even need a power supply when controlled through the computer. It is therefore the perfect solution for anybody seeking a top quality portable audiometer.

**COMPUTER CONTROL**

Connect Piccolo to the USB port of your computer and you are ready to go. There is no need to install a driver or even plug in a power supply. You can control Piccolo through the Maestro's amazingly user-friendly interface, available also for Noah software.



**ITALY**

INVENTIS S.R.L. A SOCIO UNICO  
CORSO STATI UNITI, 1/3  
35127 PADOVA  
TEL.: +39.049.8962.844

**NORTH AMERICA**

INVENTIS NORTH AMERICA INC.  
2503 S WASHINGTON AVE #586  
TITUSVILLE FL, 32780  
TOLL FREE: +1.844.683.6847

**FRANCE**

SYNAPSYS SAS  
2 RUE MARC DONADILLE  
HÔTEL TECHNOPTIC  
13013 MARSEILLE  
TÉL.: +33.4.91.11.75.75